


3880 MRI Patient Monitor System Reference Chart

SELECT PATIENT PROFILE

a. Press  Button on the 3880Q Monitor

b. Press Store/Recall Setups Button



c. Press stored setting and then press Recall Button



MONITOR BUTTONS



Monitor Settings Button
(Access main menus)



StandBy (Alarm Off) Button

(Indefinitely pauses audio and visual alarms, Terminates NIBP cycle)
(Does NOT put monitor to sleep)



Trend Button
(Access vital sign trends)



NIBP (Star/Stop) Button (Activates NIBP measurement/cycle)
Pressing button for 3 seconds starts stat reading



Record Button
(Start/Stop Base Station printer)






Alarm Silence (Audio Alarm Off) Button

(Silences audio alarms for 120 seconds, visual alarms continue to flash)

NON-INVASIVE BLOOD PRESSURE (NIBP)

a. Select NIBP cycle time




b. Press  Button to initiate NIBP cycle

* The monitor will not auto cycle NIBP when in Stand-By Mode

* Neonate or Infant size NIBP cuffs can only be used with Neonate Mode


WIRELESS POD CHANNEL CONNECTIVITY



a. Press  to power ON POD.

b. Verify Connected light is illuminated.

c. If not Connected:

- Press  to sync POD with Monitor.
- Continue pressing until **Connected light is ON.**
- ALT Channels light will illuminate if on CH 5 - 8.


MONITOR ALARM OVERVIEW



On Screen Alarms		
Alarm Priority	Physiological Alarms	Technical Alarms
HIGH	Possible Life-Threatening Situation RED Flashing Dome Light and Parameter Numerics Audible Sound.	Indicates a serious device related issue. Same Visuals/Audio as Physiological Alarms Technical Message Displayed.
MEDIUM	Possible Serious Non-Life-Threatening Situation YELLOW Flashing Dome Light & Parameter Numerics Audible Sound.	Device related issue or improper Operation Same Visuals/Audio as Physiological Alarms Technical Message Displayed.
LOW	Vital Sign Abnormality- Requires Clinical Attention BLUE Flashing Dome Light and Parameter Numerics Audible Sound.	Device related issue may compromise monitoring functions. Same Visual/Audio: Technical Message Displayed

Alarm Examples:	
High Priority MESSAGE	X
Med Priority MESSAGE	X
Low Priority MESSAGE	X

SpO₂ PULSE OXIMETRY APPLICATION

Adult & Pediatric Disposable SpO₂ Wrap

a. Apply disposable SpO₂ wrap on patient's finger or toe with  symbol over nail bed. Wrap wings of SpO₂ wrap around patient's finger or toe.



b. Insert white tipped cable with **red light** into opening on grip with  symbol.

c. Insert black tipped cable tip into grip's bottom opening with  symbol.

d. Fiber-optic tips must face each other.




Neonate and Infant/Adult Disposable SpO₂ Wraps

a. Apply disposable wrap with  symbol over nail bed, top of foot, or hand.

b. Insert white tipped cable with **red light** into opening on  grip.

c. Insert black tipped cable tip into grip's bottom opening marked by with  symbol.

d. Fiber-optic tips must face each other.



Reusable SpO₂ Grips

a. Apply grip with **number** over nail bed and facing away from patient.

b. Insert white tipped cable with **red light** into opening with **number** on nail bed.

c. Insert black tipped cable into grip's bottom opening on finger pad.

d. Fiber-optic tips must face each other.



SCAN FOR TRAINING VIDEOS

ECG



SPO2



IBP



CO2



NIBP



SETUP



ECG TRBL




PLACE WIRELESS PODS INSIDE PROTECTIVE **POD JACKETS**.





ECG ELECTRODE APPLICATION

Place **electrodes & lead-wires** on patient as follows:



Separate electrodes for Adult ECG application



Pediatric



Infant

Leadwire Management

- Route PODs & leadwires away from RF Coil. Place PODs in "Line of Sight" with Monitor.

ECG Trouble Shooting:

1. QRS waveform amplitude should be larger than the **white bracket**.
2. Ensure 3880 Monitor's MRI Filter is set to **MRI or MRI 2** to filter out RF interference on waveform.



3. Change **Leadview** if needed.
4. Prep Skin with Nuprep and apply new electrodes.

Three Lead Placement

- a. Remove **V_G** & **RL_N** wires from ECG POD
- b. Place electrodes as follows



Potential Causes of Waveform Noise

1. Weak ECG amplitude due to placement
2. Non-MR Safe or expired ECG electrodes
3. Alcohol based skin prep products
4. Damaged ECG leadwires
5. POD or ECG Leadwire under MRI coil
6. New gradient sequence

INVASIVE BLOOD PRESSURE

a. Power ON and Connect **transducer** cable to P1 or P2 on ipPOD. **DO NOT ATTACH TRANSDUCER TO PATIENT SKIN.**

Transducer must be maintained at applicable level for invasive pressure per hospital policy.

*Must use MR compatible transducer cable (not sold by IRadimed) with 3883-1 Universal ipPOD.



3883-1(Universal)*



3883-2 (Truwave)




3883-3(Transpac)

b. Keep transducer and ipPOD **OUTSIDE** of MRI bore.

c. Select **IBP** section on monitor's screen. Select IBP label and set **Scale** as applicable.

Turn transducer stopcock "OFF TO PATIENT/OPEN TO AIR."

d. Remove non-vented cap from stopcock and Press "**CAL ZERO**" in IBP section of monitor screen corresponding with **IBP1** or **IBP2** as applicable.



e. When **00/00** is displayed in the IBP parameter section, turn stopcock "OPEN TO PATIENT/ CLOSED TO AIR."

Replace non-vented cap on stopcock.

EtCO₂ FROM MONITOR (INTERNAL)

a. Attach EtCO₂ Sample Line with Disc Filter to CO₂ port on monitor

b. Connect Sample Line to patient's airway adapter, nasal cannula, or mask



Factors that may affect accuracy:

- Leaks, kinks, or moisture in Sample Line or Filter
- Airway pressure
- Use of masks or extra extension tubing.
- Incompatible accessories for breathing style



1841 micron disc filter is single patient use only.

EtCO₂ FROM 3886 MULTI-GAS MODULE (EXTERNAL)




MR Conditional up to 600 Gauss

Multi-Gas Module provides CO₂ waveform and anesthetic agents measurements

a. Power Multi-Gas Module **ON** using switch on rear side

b. Insert Nomoline into Multi-Gas Module



Ready to Measure Gases



c. Attach Sample Line between Nomoline & patient's airway adapter/cannula

CO₂ OCCLUSION MESSAGE: Check and/or replace Sample Line or Nomoline.

1850 nomoline can be reused for up to 14 days.